

Thermo-Cell-Transporter 3018

Qualität – Made in Germany

Zum sicheren Transport

For safe transport

- konstant 37°C
- zuverlässige Thermostatisierung
- netzunabhängig
- vielseitig einsetzbar

- constantly 37°C
- reliable temperature regulation
- independent of mains connection
- versatile use



Technische Daten

Gehäuse:

- Abmessungen 111 x 111 x 215
(B x T x H in mm)
- Gewicht 2,4 kg

Temperatur:

- 37°C konstant
(werkseitig fest justiert; Sonderanfertigung auf Wunsch im Bereich von 32°C bis 42°C möglich)

Regelgenauigkeit:

- $\pm 0,3^{\circ}\text{C}$

Betriebszeit:

- Netzteil \Rightarrow Dauerbetrieb
- Kfz-Spannung \Rightarrow Dauerbetrieb
- Akku-Betrieb \Rightarrow ca. 5 Std.*

Aufheizzeit/Ladezeit:

- 22°C – 37°C \Rightarrow ca. 1 Std.*
- Akku-Ladezeit \Rightarrow ca. 3 Std.*

Anschluss:

- an Buchse „+12V Ext.“ \Rightarrow 11V...15V

Lieferumfang:

- 1 Grundgerät
- 1 Netzteil
- 1 Kfz-Kabel

Zur Verwendung mit:

- 8 Falconröhrchen (Rundboden, 13 ml) oder
- 1 Multischale (66 mm x 66 mm) oder
- 1 Gewebekulturschale (\varnothing 60 mm) oder
- 1 Petrischale (\varnothing 60 mm)

* bei einer Umgebungstemperatur von 22°C



Technische Änderungen vorbehalten

Technical Data

Housing:

- Dimensions 111 x 111 x 215
(w x d x h in mm)
- Weight 2.4 kg

Temperature:

- 37°C constant
(factory-aligned; other specifications available upon request in the range of 32°C-42°C)

Accuracy of regulation:

- $\pm 0.3^{\circ}\text{C}$

Operation time:

- Power supply \Rightarrow continuous operation
- 12V car voltage \Rightarrow continuous operation
- operation with rechargeable batteries \Rightarrow approx. 5h*

Heating/Charging time:

- 22°C – 37°C \Rightarrow approx. 1h*
- charging time of batteries \Rightarrow approx. 3h*

Connection:

- at socket „+12V Ext.“ \Rightarrow 11V ...15V

Delivery volume:

- 1 Basic unit
- 1 Power supply
- 1 vehicle connection cable

For use with:

- 8 Falcon tubes (round bottom, 13 ml) or
- 1 Four-well dish (66 mm x 66 mm) or
- 1 Tissue culture dish (\varnothing 60 mm) or
- 1 Petri dish (\varnothing 60 mm)

* at 22°C ambient temperature



Specifications subject to change without notice

ÄM72_08_Rev. 0