

## Accessories for active sterile humidity supply of the CO<sub>2</sub> incubators C60, C200 & C201

### Advantages:

- Tested water
- Increased safety with disposables
- No extensive cleaning and sterilization of the external water reservoir
- Easy handling

**AEROp<sup>®</sup>part HUM**  
sterile and  
pyrogen-free water  
1000 ml

REF 16150

**Intrafix<sup>®</sup> Air (B. Braun)**  
I.V. administration set  
tubing length 150 cm

REF 15417

**Magnetic hook**  
diameter 48 mm

REF 15595



Rev. 6\_02/2020