

## Vibration Technology

Seismion GmbH Berliner Str. 8 30457 Hannover Germany

sales@seismion.com www.seismion.com

## Reactio

## Feedback Control Vibration Isolator

Made in Germany



## **Key Features**

- High performance, low price and different sizes offer cost efficient access to active vibration isolation for your applications
- Actively stabilized top-plate in 6-degrees of freedom, counter forces disturbances from applications, i.e. scanning X/Y - stages
- Super flat design, no configuration required and easy to use
- Meets highest vibration criteria: VC-F @ 1 Hz and VC-G @ 2 Hz
- Robust electronic feedback control handles even strong vibration and disturbances without interrupting or stopping isolation



