

Cryo Watch

Quality – Made in Germany

Temperature guard for monitoring of temperature inside LN₂ cooled cryo storage tanks



Features:

- Temperature measuring up to - 196 °C
- Temperature is displayed continuously
- Temperature sensor inside the cryo tank
- Acoustic, optical and remote alarm if set point is breached
- Programmable alarm set point from - 196 °C to + 40 °C
- Remote alarm active in case of power failure
- Usable as universal temperature monitoring system

Technical data:

- Temperature range: - 196 °C to + 40 °C
- Accuracy: +/- 2 °C
- Dimensions: 200 x 90 x 100 mm
- Power supply: 230 V AC, 50/60 Hz, 5 W
- Safety class II



REF 13852



Rev. 5_01/2017