

Fermacidal D2®



- Without alcohol
- Fast acting
- MEA and HSSA tested

Product:

Spray bottle 50 ml	REF	13045
Spray bottle 250 ml	REF	13048
Bottle (screw cap) 500 ml	REF	15213
Spray bottle 1 l	REF	13044
Canister 5 l	REF	13043
Canister 10 l	REF	15101
Disinfectant wipes, 120 pieces	REF	13046

Advantages:

- Ready-to-use disinfectant for surfaces and instruments
- Cleaning effect
- Long-term effective
- Smell-removal and neutralisation of odour-causing bacteria
- Also suitable for disinfection of CO₂ incubators

Microbiological efficacy:

- Bactericidal, e.g.: *Salmonella*, *Mycobacterium tuberculosis*
- Fungicidal, e.g.: *Trichophyton mentagrophytes*
- Selective virucidal, e.g.: *Hepatitis B*, *HIV*, *Rotavirus*, *Influenza A*, *Corona viruses*

(as specified by the manufacturer)

Features:

- Quaternary ammonium compounds
- Without alcohol and aldehyde
- Odourless
- Non-volatile
- Does not irritate the skin
- Does not stain
- Dissolves dried blood stains
- Good compatibility of materials like metal, rubber and plastic
- VAH listed and BauA registered

(as specified by the manufacturer)



Rev. 6_06/2020