

## Catheter for hysterosalpingo sonography (HSO catheter set)

Quality – Made in Germany

**Sterile disposable balloon catheter for intrauterine hysterosalpingo contrast sonography**

**Available in two packaging units:**

• 1 piece: PZN 9116253    REF 14898

**Consisting of:**

1. HSO catheter
  - Side opening for contrast medium
  - Fixation and occlusion balloon
  - Disposable injection nozzle for contrast medium
  - Stylet
  - Stopcock
2. Filling syringe



Contains natural latex

Rev. 4\_05/2018