

Hot Plate 100

Quality – Made in Germany



Hot Plate for laboratory and scientific use

Versatile use between 27 °C and 100 °C: Heating, drying, stretching

Features:

- Temperature control between 27 °C and 100 °C
- Homogeneous temperature distribution
- Planar surface for optimal heat transfer to heating blocks, culture dishes and slides
- Short heating up time
- Compact, flat, with carry handles
- Space saving due to separate control unit

Technical data:

- Dimensions: 314 x 214 x 23 mm (w x d x h)
- Weight: 2.9 kg
- Heating up time (approx.):
37 °C: 2 min., 100 °C: 9 min.
- Power supply:
100 – 240 V AC, 50/60 Hz, 140 W
- Temperature adjustable between 27 °C and 100 °C
- Safety class I

Scope of supply:

- Hot plate
- Controller with mains adapter

REF 15049



Rev. 4_02/2020