

## NEW: Hot Plate A4

Quality – Made in Germany



### Hot Plate for laboratory and scientific use

#### Features:

- Robust and planar glass surface
- Homogeneous temperature distribution
- OLED display for intuitive and easy-to-use user menu
- Functional buttons with haptic feedback via vibration (disengageable)
- Extremely flat and space saving
- Cleaning and disinfection with all common detergents

#### Technical data:

- Dimensions: 300 x 210 x 15 mm (w x d x h)
- Usable surface area: 255 x 140 mm
- Temperature stability:  $\pm 0.2$  °C
- Temperature homogeneity:  $\pm 0.3$  °C
- Weight: 1050 g (w/o mains adapter)
- Heating up time: approx. 10 min
- Adjustable temperature range: 27-45 °C (2 °C above ambient temperature)
- Power supply: 100 – 240 V AC, 60 W
- Cable length: Hot Plate to mains adapter 1.8 m, mains adapter to plug 1.8 m

(Specifications subject to change without notice.)

#### Scope of supply:

REF 16481 Hot Plate A4 with mains adapter and cable

**CE** Rev. 2\_02/2018