

Hot Plate A4

Quality – Made in Germany



Hot Plate for laboratory and scientific use

Features:

- Robust and planar glass surface
- Homogeneous temperature distribution
- OLED display for intuitive and easy-to-use user menu
- Functional buttons with haptic feedback via vibration (disengageable)
- Integrated Power Switch
- Extremely flat and space saving
- Cleaning and disinfection with all common detergents

Technical data:


- Dimensions: 300 x 210 x 15 mm (w x d x h)
- Usable surface area: 255 x 140 mm
- Temperature stability: ± 0.2 °C
- Temperature homogeneity: ± 0.3 °C
- Weight: 1050 g (w/o mains adapter)
- Heating up time: approx. 10 min
- Adjustable temperature range: 27 – 45 °C (starting at 2 °C above ambient temperature)
- Power supply: 100 – 240 V AC, 60 W
- Cable length:
 - Hot Plate to connection point on power supply unit 0.3 m
 - Power supply unit to plug 1.8 m

Scope of supply:

- Hot Plate with mains cable REF 16481

Additional equipment:

- Retractor frame REF 17958



(Specifications subject to change without notice.)