

Incubator-Filter-Box

Quality – Made in Germany



Separate unit for filtering the atmosphere in incubators

Features:

- Suitable for all cell culture incubators with minimum 60 l volume
- Eliminates particles, VOCs and other impurities in a very short time
- External controller outside the incubator
- No negative influence on the atmosphere inside the incubator
- Flexible use

Optional filters:

Gas filter A1 14548

Composite filter A2-P3D 15415

Technical data:

- Composite filter is preinstalled = standard filter (to be replaced every 3 months)
- Operation under regular incubator atmosphere
- Dimensions: diameter 10 cm / height 12 cm
- Weight: approx. 600 g (when using composite filter)
- Switches on automatically for 5 - 6 minutes; starts again after a break of 40 minutes

14526



Rev. 7_04/2017