

Icepur®



- **Disinfectant cleaner concentrate**
- **Free of alcohol**
- **Free of aldehydes**

Product:

Grip bottle 2 l

REF 15422

Canister 10 l

REF 15421

Advantages:

- For combined wipe disinfection and cleaning of surfaces and inventory
- For mask disinfection and cleaning in case of epidemics as well
- Free of any toxic or volatile components
- Protein and fat dissolving, neutralizes odour-producing bacteria

Microbiological efficiency:

- Bactericidal,
e.g.: *Salmonella*, *Mycobacterium tuberculosis*
- Fungicidal,
e.g.: *Trichophyton mentagrophytes*
- Selectively virucidal,
e.g.: Hepatitis B, HIV, Rotavirus

(as specified by the manufacturer)

Features:

- Quaternary ammonium compounds
- Without alcohol and aldehydes
- Odourless, neutralizes odour-producing bacteria
- Non-volatile
- Cleans and disinfects in one step
- Dissolves dried blood stains
- Leaves no stains
- Protein and fat dissolving
- Good compatibility with metals, rubber and plastics
- BauA registered

(as specified by the manufacturer)



Rev. 0_05/2020