


LaboDat+

Quality – Made in Germany



Software for documentation of incubation processes

Figure: Main window LaboDat+

Features:

- Display of actual and stored values
- Protected record of all measurements and alarms
- Visualization in diagrams or numerical values
- Compatible to Labotect devices with RS485 interface*
- Data export possible
- Language selection English / German
- Compatible up to Windows 10

Hardware:

- USB interface converter

REF 15125

*depending on software version of the main board

Advantages:

- Documented record of incubation conditions
- Connection for up to 16 devices
- USB interface
- Report generation
- Data storage on server possible

Compatible devices:

- ➔ C16
- ➔ C60
- ➔ C200
- ➔ Labo C201
- ➔ Labo C-Top
- ➔ Labo Gas Mixer
- ➔ Gasmonitor




Figure: LaboDat+ evaluation of single device