

Retractor frame

Quality – Made in Germany



Picture:
Retractor frame with Hot Plate 062

Accessory for:

- HotPlate 062
- HotPlate A3

- For variable fixing of objects to magnetic holders
- Stainless steel (magnetic)

Retractor frame for HotPlate 062

Inner dimensions: 21.0 x 15.0 cm
Outer dimensions: 27.7 x 21.7 cm

REF 15943

Retractor frame for HotPlate A3

Inner dimensions: 34.0 x 22.0 cm
Outer dimensions: 42.4 x 30.4 cm

REF 15942

Scope of supply:

- Retractor frame
- Silicone base
(Magnetic holders not in scope of supply)



Rev. 1_01/2017