

## Retractor frame

Quality – Made in Germany



Picture:  
Retractor frame with Hot Plate A3



Picture:  
Retractor frame with Hot Plate A4

### Accessory for:

- Hot Plate 062
- Hot Plate A3
- Hot Plate A4

- For variable fixing of objects to magnetic holders
- Stainless steel (magnetic)

### Scope of supply:

- for Hot Plate A3:
  - Retractor frame
  - Silicone base
- for Hot Plate A4:
  - Retractor frame with integrated silicone feet

(Magnetic holders not in scope of supply)

### Retractor frame for Hot Plate 062

Inner dimensions: 21.0 x 15.0 cm  
Outer dimensions: 27.7 x 21.7 cm

REF 15943

### Retractor frame for Hot Plate A3

Inner dimensions: 34.0 x 22.0 cm  
Outer dimensions: 42.4 x 30.4 cm

REF 15942

### Retractor frame for Hot Plate A4

Inner dimensions: 24.7 x 13.2 cm  
Outer dimensions: 30.4 x 21.4 cm

REF 17958



Rev. 3\_05/2021