

## Thermo Cell Transporter 3018

Quality – Made in Germany



Transport incubator for safe and versatile use

### Features:

- Factory-aligned temperature at 37 °C
- Different temperature possible (to be advised before order)
- Continuous operation with power supply, car voltage or integrated accumulator during transport
- Operation with batteries up to 5 h\*
- Compact design

\* at 22 °C ambient temperature

### Technical data:

- Dimensions: 111 x 111 x 215 (w x d x h in mm)
- Weight: 2.4 kg
- Power supply: 100 – 240 V AC, 50/60 Hz, 60 W  
12 V DC, 1.5 A car connection cable
- Temperature: between 32 °C and 42 °C available (to be advised before order)
- Accuracy of regulation:  $\pm 0.3$  °C
- Operation time:
  - with power supply or 12 V car voltage continuously
  - with accumulator approx. 5 h\*
- Heating up time 22 °C to 37 °C: approx. 1 h\*
- Charging time of accumulator: approx. 3 h
- For use with 8 x 13 ml Falcon tubes (REF 13881) or 2 x 13 ml and 8 x 5 ml Falcon tubes (REF 16518)
- Safety class II

### Scope of supply:

- Basic unit
- Power supply unit
- Car connection cable

REF 13881 for 8 x 13 ml Falcon tubes

REF 16518 for 2 x 13 ml and 8 x 5 ml Falcon tubes



Rev. 5\_02/2020